











Endoscopic forceps

Laparoscopic instrument

BICARmed technology guarantees maximum effectiveness on all instruments and dirt types, even the most adherent contaminant will be removed.

BICARmed TECHNOLOGY IS DISTRIBUTED IN FRANCE BY



PRE-CLEANING OF RMD WITH **BICARBONATE**

INNOVATIVE TECHNOLOGY

> Which guarantees you high efficiency, before going through the washer-disinfector. Thanks to the precise projection of a solution of sodium bicarbonate, compressed air and water at low pressure, without corrosion or abrasion.

In a vacuum cabin, the protected operator loads the MD to be cleaned through the side door. The washing nozzle strikes the surface of the instrument evenly and continuously (98.2% efficient).

- > Possible installation of the machine in Sterilization or in the operating room.
- > Provision of a machine at CH Rennes